



Infant Patient ID Solutions



Most Labor & Delivery and NICU traditionally use a serial numbered mother/father/baby wristband set, a method which provides important benefits such as matching the baby to its parents and the ability to identify the baby immediately after delivery. Often, handwritten information has also been used. However, today's healthcare systems require patients to have an attached barcode identifier for a variety of applications such as medicinal administration, lab work, sleeping and feeding schedules, etc. A common concern among hospital staff and infant caregivers is a proper way to quickly and efficiently "barcode a baby".

Specific Caregiver Requirements

Hospital caregivers have some very specific identification tag requirements for these special patients:

- Must be secure and not removable.
- Must enable caregivers to get an easy scan of the barcode without jostling the baby.
- Must be easy to generate and replace at the department level.
- Must be able to reuse if serial numbered band is removed for a procedure.



Our Newborn Identification Tags

There is no need to replace your current method. You can easily add our Newborn Identification Tags to the serial numbered bands you use today.

- Are secure and cannot be removed and reused.
- Hang flat so they are easy to scan.
- Are printable on both sides for easy scan capture.
- Are durable, yet comfortable.
- Can be printed on an inexpensive "plug and play" thermal printer.
- No need for your IT department to write a custom interface.



**Printable on
Both Sides**




www.IdentificationSystemsGroup.com

ISG™ Newborn Identification Tags

ISG™ ContourCare® Newborn ID Tags

Achieve positive identification for newborn patients using ISG™ ContourCare® Newborn Identification Tags. The unique, double-sided design and adhesive-free strap provide a flat barcode surface for easy access and scanning without disturbing the patient. The tags are made from soft material with rounded corners so they are gentle on your most sensitive patient's skin.



- Double-sided tag with strap allows scanning at any position, without disturbing the patient
- Easy-to-use and attaches to any wristband
- Durable, lightweight and comfortable design
- Latex-free and phthalate-free
- Supports text, barcodes and photos
- Compatible with most direct thermal printers

Product Number	Core Size	Wound	Tags per Box	Information Area
8525-83-ISM	1"	In	500	1.6875" x 1.84375"
8526-83-ISM	1"	Out	500	1.6875" x 1.84375"
8527-83-ISM	1 - 1/2"	In	500	1.6875" x 1.84375"
8528-83-ISM	1 - 1/2"	Out	500	1.6875" x 1.84375"

For more information about ContourCare® healthcare patient identification solutions from the ISG, or related printer models and supplies, please visit our website at www.IdentificationSolutionsGroup.com or call us at 888-964-6482 to speak to your local ISG Dealer.



Identification Systems Group

— The ISG™ is a nationwide network of local experts in identification, security, and card personalization, providing high quality, cost-effective solutions backed by local support and the strength of our Professional Services Certification program. Each company works together to provide seamless nationwide support and collaboration in the identification and card issuance industry.

www.IdentificationSystemsGroup.com