



Key Features

- Double-sided tag with strap allows scanning at any position
- Easy-to-use & attaches to any wristband
- Durable, lightweight & comfortable design
- Latex-free & phthalate-free
- Supports text, barcodes & photos
- Printable on both sides
- Compatible with most direct thermal printers

ContourCare® Newborn Identification Tags

Infant patient ID solutions
that meet specific caregiver
requirements

Most Labor & Delivery and NICU traditionally use a serial numbered mother/father/baby wristband set, a method which provides important benefits such as matching the baby to its parents and the ability to identify the baby immediately after delivery. Often, handwritten information has also been used. However, today's healthcare systems require patients to have an attached barcode identifier for a variety of applications such as medicinal administration, lab work, sleeping and feeding schedules, etc. A common concern among hospital staff and infant caregivers is a proper way to quickly and efficiently "barcode a baby".

Hospital caregivers have some very specific identification tag requirements for these special patients:

- Must be secure and not removable
- Must enable caregivers to get an easy scan of the barcode without jostling the baby
- Must be easy to generate and replace at the department level
- Must be able to reuse if serial numbered band is removed for a procedure

For more information about newborn healthcare patient identification solutions, related thermal printer models and supplies, contact the ISG or your local dealer.




ContourCare® Newborn ID Tag Versions

Product Number	Core Size	Wound	Tags per Box	Information Area
8525-83-ISM	1"	In	500	1.6875" x 1.84375"
8526-83-ISM	1"	Out	500	1.6875" x 1.84375"
8527-83-ISM	1-1/2"	In	500	1.6875" x 1.84375"
8528-83-ISM	1-1/2"	Out	500	1.6875" x 1.84375"

Achieve positive identification for newborn patients using ContourCare® Newborn Identification Tags. The unique, double-sided design and adhesive-free strap provide a flat barcode surface for easy access and scanning without disturbing the patient. The tags are made from soft material with rounded corners so they are gentle on your most sensitive patient's skin.

There is no need to replace your current method. You can easily add our Newborn Identification Tags to the serial numbered bands you use today.

- Are secure and cannot be removed and reused.
- Hang flat so they are easy to scan.
- Are printable on both sides for easy scan capture.
- Are durable, yet comfortable.
- Can be printed on an inexpensive "plug and play" thermal printer.
- No need for your IT department to write a custom interface.



About the ISG

The Identification Systems Group (ISG) is a nationwide network of local experts in identification, security, tracking and card personalization technologies, providing high quality, cost-effective solutions backed by local support and the strength of our Professional Services Certification program. Each member company works together to provide seamless support and collaboration in the identification and issuance industries across the USA and Canada.

www.IdentificationSystemsGroup.com

